


## AZDye 405 Alkyne

<b>Product Name:</b>	AZDye 405 Alkyne
<b>Product No.:</b>	1067
<b>Molecular Weight:</b>	553.53 protonated 857.11 triethylammonium salt
<b>CAS:</b>	n/a
<b>Solubility:</b>	Water, DMSO, DMF, MeOH
<b>Appearance:</b>	Grey solid
<b>Shelf life:</b>	At least 12 months at -20°C
<b>Shipping:</b>	At ambient temperature
<b>Storage:</b>	Upon receipt store at -20°C

### Chemical Structure



### Optical Properties

$\lambda_{\text{abs}}$	402 nm
$\lambda_{\text{abs}}$	424 nm
$\epsilon_{\text{max}}$	73,000 M <sup>-1</sup> cm <sup>-1</sup>
$\Phi_{\text{fl}}$	n/a
CF <sub>260</sub>	0.78
CF <sub>280</sub>	0.70
Spectrally Similar	Alexa Fluor®405
Dyes	CF® Dye 405, DyLight® 405

DyLight is a trademark of Thermo Fisher Scientific, Inc. Alexa Fluor is a trademark of Thermo Fisher Scientific, Inc. CF Dyes is a trademark of Biotium Inc.

### Absorption and Emission Spectra

